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Elekta Synergy[®]

Digital accelerator for advanced IGRT

*Setting the standard
for confident care*


Elekta Synergy[®]

Inspiring clinical confidence

Elekta Synergy was the first advanced digital linear accelerator to introduce soft tissue volumetric Image Guided Radiation Therapy. With a suite of advanced guidance tools, imaging of soft tissue at the time of

treatment allows efficient verification of tumor and critical structure position, providing increased confidence in dose placement.


Elekta Synergy[®] Platform

Foundation for the future

As the foundation of the Elekta Synergy family of digital accelerators, Elekta Synergy Platform has been designed to grow with your practice. It introduces real-time patient position verification, to ensure confidence in patient set-up accuracy, and efficient guidance of dose placement.


Through the addition of kV imaging to its current IGRT technology, the Synergy Platform can be easily upgraded to Elekta Synergy, allowing centers to offer a choice of imaging modalities best suited for the patient's need.

Patient-centric

The application of evidence-based medicine allows informed decision making on treatment techniques, outcomes and performance. From patient-centric treatment planning solutions and MOSAIQ®

comprehensive patient management information system, Elekta Synergy® gives users access to the entire Oncology chart in a single integrated solution.

Oncology Information Systems

At the heart of Oncology Management, MOSAIQ® image-enabled electronic medical record provides tools to streamline the radiotherapy process as a whole, while delivering secure access to the patient information and images that drive clinical decision-making.

Treatment Planning

Elekta's patient-centric planning solutions provide a broad range of leading edge applications and comprehensive workflow solutions, supporting the requirement for increasingly tighter coupling between treatment planning, delivery workflow management and review.

Unmatched Clinical Confidence

With IGRT becoming an integral part of the standard of care for many patients, Elekta Synergy® not only provides the ability to deliver tailored imaging through a variety of imaging modalities, but includes many unique features,

such as the industry's largest IGRT field of view, that enable you to deliver personalized patient care throughout the entire radiotherapy treatment process.

Imaging

For unmatched clinical confidence, Elekta Synergy delivers a choice of 2D MV, 2D, 3D and 4D kV image guidance for visualization and tumor verification. With a suite of advanced image guidance tools, Elekta Synergy provides clinicians with the ability to provide tailored image guidance to meet the needs of each individual patient.


Image Guided Treatment Management

Through coordinated image guidance and treatment delivery Elekta Synergy® provides a synchronized workflow for Radiation Oncology. Fundamentally image-enabled, MOSAIQ provides the flexible treatment management solution designed to support the complexities of image-guided radiation therapy.

Digitally Controlled, Clinically Focused

Accuracy and reliability of treatment delivery performance are imperative to be able to deliver individual patient care with confidence. Powered by Elekta's 6th generation of integrated digital control, Elekta Synergy[®] seamlessly combines advanced

guidance though IGRT with flexible treatment delivery across a wide number of techniques. With the market leading isocentric clearance of 45cm, Elekta Synergy provides significantly more freedom to deliver non-coplanar therapies.

Digital Control

Through simple, intelligent coordination of multiple linac functions, Elekta Synergy is powered by digital control that orchestrates the linac's myriad of operations to ensure that any treatment is delivered safely, accurately, quickly and intelligently.


Treatment Delivery

Individualized patient care is possible through a range of techniques from 3D conformal techniques, through static and dynamic IMRT to VMAT. With these advanced delivery techniques, you no longer have to choose between speed and accuracy.

Flexible Beam Shaping

Supporting the range of treatment delivery techniques, Elekta Synergy supports a range of integrated multileaf collimators, allowing users to choose the beam shaping that meets the need of their clinic. Through innovative

design, Elekta Synergy provides the ability to reduce dose to critical structures while meeting today's clinical desire to deliver higher doses to the target.


Beam Shaping

With a range of advanced integrated multileaf collimators, Elekta provides excellent beam-shaping capabilities across the range of delivery techniques. Through fully integrated digital control with continuous real-time optical verification, accurate placement of all leaves allows for faster, safer and more accurate delivery.

Active Motion Management

Delivery of high doses associated with today's radiation therapy techniques requires not only accurate reproduction of patient positioning, but also management of internal patient motion during both

imaging and treatment delivery. Elekta's clinically efficient solutions address the challenges of patient immobilization and actively manage motion, facilitating margin reduction and dose escalation.

Motion Management and Patient Positioning

Through a versatile range of positioning and immobilization solutions, optimal treatment setup is achieved through innovative designs which focus on the patient experience as a whole, combining patient comfort, safety and reproducibility.


Respiratory Motion Management

Addressing target localization associated with respiratory motion management can be challenging. Elekta Synergy's 4D anatomically correlated motion management capabilities allow clinicians to visualize respiratory motion and account for baseline shifts.

Expanding Clinical Practice

More and more clinics are undertaking practice expansion in stereotactic treatment techniques. Supporting you in this area, Elekta Synergy has been designed to be compatible with the Elekta range of stereotactic

accessories, providing you with the option to introduce stereotactic practice into your clinic simply, efficiently and cost-effectively.

Stereotactic Accessories


With a range of stereotactic accessories, including add-on micro-MLCs, precise, small-field delivery enables implementation of stereotactic radiosurgery and stereotactic radiation therapy.

Quality Assurance

Completing the picture; Elekta offers a full range of QA tools for fast accurate verification of performance, providing the confidence that treatments are delivered as planned.

XVI Intra-fraction Imaging

Monitor and manage internal motion

*2D, 3D & 4D soft tissue imaging
during treatment delivery*






Building on proven excellence in advanced image guidance

Elekta was the first to develop 2D, 3D and 4D kV soft tissue imaging that allowed visualization of not only the tumor but also the surrounding critical structures at the time of treatment. In recent years Elekta has taken this further by moving from soft tissue imaging to true treatment guidance. Today with XVI intra-fraction imaging, Elekta takes another leap forward in IGRT towards greater excellence and confidence in patient care by providing the tools to monitor and manage internal motion.

XVI intra-fraction imaging is a part of XVI 5.0 and builds on the proven ability and excellence of XVI from Elekta. With the ability to image during treatment delivery, it allows for opportunity to reduce treatment time-slots to maximize clinical efficiency.

Combining imaging and treatment delivery also reduces the likelihood of patient movement and changes in internal organ position during the treatment session. This means that patient care is further enhanced while allowing flexibility to clinicians to provide a workflow that is specific to each patient.

Why XVI Intra-fraction Imaging?

2D Intra-fraction Imaging

This new tool allows visualization of soft tissue images during treatment either as static or fluoroscopic-like images. This real-time monitoring provides opportunity to monitor and correct for motion during treatment. With this 2D intra-fraction imaging approach critical structures can be protected by utilizing real-time images that can be interpreted using the annotation overlay tool. This functionality allows structures drawn onto the reference image to be overlaid on the real-time image that is being captured during treatment.

3D Intra-fraction Imaging

The ability to acquire 3D VolumeView™ images during rotational treatment delivery including VMAT, provides a post treatment record of patient position verification. This allows for online and offline correctional strategies to be implemented utilizing information from the exact moment of treatment. With this additional information, confidence can be further increased by accounting for changes before the next treatment session or even between treatment fields.

4D Intra-fraction Imaging

Symmetry™ introduced a new approach in the management of respiratory motion. With the new XVI and its ability to take images during delivery, this approach can provide a verification of real tumor motion as treatment delivery is being undertaken. With moving tumors being a complex challenge in radiation treatment, this confirmation of real-time tumor movement during delivery provides added confidence especially when looking to escalate doses.

Workflow Reinvented

Utilizing the power of MOSAIQ® image-enabled oncology information system (OIS), XVI 5.0 also provides greater flexibility in the use of image-guided workflows to enable distributed workflows for clinical efficiency and optimized use of resources. The advanced integration with MOSAIQ means that information regarding each patient's imaging scan, setup and correction are available as an integral part of the patient's electronic medical record.

Intra-fraction 2D, 3D & 4D for management of internal motion

Verify anatomical position and motion during treatment

Increased confidence for dose escalation

Distributed workflows with MOSAIQ for further clinical efficiency

XVI 5.0 is not available in all regions. Please contact your local representative for more details.

Agility

ELEKTA

High-resolution full field MLC

Leaf motors drive each of the 160 leaves to the position determined by the integrated digital control system

Dynamic leaf guides ensure seamless field delivery across full 40cm field width and contribute to industry leading leaf speed

Rubicon provides real-time optical leaf positioning

160 leaves, 5mm wide at the isocenter, providing beam shaping over the full 40 x 40cm field size.

Rubies attached to each individual leaf tip. They fluoresce in the near-UV light and emit light in the near-IR wavelength.

Video camera picks up light in the near-IR range. The image created is used by the digital control system to monitor leaf position and ensure positional accuracy.

Near-UV source floods the head with light in the near-UV spectrum.

Sculpted diaphragms provide shielding in the axis perpendicular to the leaf direction. Leaves position themselves intelligently behind the Y shaped portion of the diaphragm.

Monaco[®] 5

Treatment planning system

Accuracy accelerated


Monaco is a next-generation treatment planning system. A comprehensive range of techniques is now supported in Monaco, including 3D conformal radiation therapy, IMRT, VMAT, stereotactic MLC and cone, 4D and Siemens mARC (rIMRT).


IMRT

Monaco includes intuitive and leading edge features for generating and evaluating IMRT plans. Biological cost functions and constrained optimization lead to predictable results with evaluation tools for where to improve the overall plan.

- Multiple cost function support
 - Biological: EUD for tumor volumes and serial organs, isoeffective volume for parallel organs
 - Physical: Dose volume based constraints for tumor volumes and organs at risk
- Multi-criterial when selected will automatically try to achieve better normal tissue sparing by tightening constraints during the optimization
- Enhanced voxel definition tools controls and patient model improve the planning workflow. Voxel definition tools allow the user to see where and how the patient anatomy is controlled by plan constraints
- Smart Sequencing[®] - Constrained optimization of deliverable segments lead to efficient step-and-shoot fields with fewer segments and a minimum of overall monitor units
 - Fluence (monitor units) and shape of segments are optimized subject to the initial prescription
 - No final dose calculation is necessary
- Segment Shape Optimization (SSO) improves plan quality and delivery efficiency through its proprietary approach to smoothing and clustering segments, then optimizing beam weights and shapes.
- The Monaco treatment planning system offers full support FFF by simply creating a machine model based on measurements from a FFF linac.

- Constrained optimization optimizes dose to targets while respecting hard constraints on organs-at-risk
- Constrained Optimization ensures that all organ at risk goals are achieved while optimizing the target objectives at the same time.
- New visual Sensitivity Analysis guides the user to quickly and easily solve the conflicts between target goals and dose-limiting constraints


- Robust algorithms
 - XVMC – The Monte Carlo simulation algorithm provides highly accurate dose calculation and is employed within the optimization loop
 - Collapsed cone beam for Monaco with GPU accelerated calculations for 3D
 - An enhanced pencil beam algorithm facilitates fast results for evaluating optimization parameters

VMAT

VMAT (Volume Modulated Arc Therapy) functionality within Monaco utilizes Segment Shape Optimization™ (SSO) which improves plan quality and decreases delivery time. Monaco can optimize single arcs with multiple rotations or multiple non-coplanar arcs simultaneously, fully supporting both convention and FFF delivery techniques, providing the flexibility and control needed for more complex treatment plans. Arc plans can be delivered with a single button push at the linear accelerator console. Gantry directions are automatically sorted, and all control points are seamlessly integrated into a single deliverable arc sequence. Monaco offers the accurate XVMC Monte Carlo dose engine for a continuous arc calculation as a single beam, instead of being limited to dose approximations that occur with a large number of discrete (control point) gantry angle positions. Monaco supports comparison of segmental IMRT with VMAT plans to choose the most appropriate for every patient.

mARC

Modulating Arc Therapy support for Siemens Linacs. Uses the same constraint based biological optimization as IMRT and VMAT and XVMC Algorithm

3D Conformal Planning

Using a graphical processor unit (GPU) for the collapsed-cone algorithm, Monaco increases the performance of 3D dose calculations even further. Monaco supports a range of modalities and techniques allowing clinicians to expand their treatment planning capabilities, including support of wedges, bolus and VMC++ Electron Monte Carlo

Stereotactic Planning

Full MLC-based planning including Apex™ micro-MLC support. Stereotactic cone support for static and arc fields using XVMC algorithm.

4D image Support

Monaco provides a comprehensive suite of capabilities that support the visualization and utilization of 4D image sets to create structures and review treatment options

Connectivity

Perform contouring, fusion, and plan review at any station or physician's workspace. Enhanced efficiency using template driven planning and report/batch printing generation.

Simulation

The integrated simulation application seamlessly performs contouring, beam placement, block/port definitions, and real-time DRRs. Multiple customizable views ensure efficient access to all required images and planning perspectives.

Contouring

The comprehensive array of drawing and editing tools from Elekta allow quick creation and expansion of structures. Includes edge detection/structure avoidance paintbrush tool and EZ-Sketch 3D auto-segmentation tool and PET SUV contouring.

Auto Fusion

Fully integrated fusion accommodates alignment of CT with CT/MR/PET/NM images for CT simulation or treatment planning utilization.

Plan Review

An intrinsic application offering a vast array of review tools and rapid off-line access for plan analysis and approval.

Rooms Eye View (REV)

Show the patient, gantry, couch, lighting, and decorations to simulate the treatment and verify patient positioning.

DICOM Integration

Integrated data transfer into MOSAIQ and other R&V systems including: Course ID, Plan Intent, Rx Site, Rx Dose, Dose Rate, Tolerance Table and Field ID.

Monaco 5.0 is not available for sale or distribution in all markets. Please contact your Elekta representative for details.

3D conformal radiation therapy, stereotactic cone, GPU-based Collapsed Cone calculation and dose service, are pending FDA 510(k) clearance.

www.elekta.com

Human Care Makes the Future Possible

Corporate Head Office:

Elekta AB (publ)
Box 7593, SE-103 93 Stockholm, Sweden
Tel +46 8 587 254 00
Fax +46 8 587 255 00
info@elekta.com

Regional Sales, Marketing and Service:

North America

Tel +1 770 300 9725
Fax +1 770 448 6338
info.america@elekta.com

Europe, Middle East, Africa, Eastern Europe, Latin America

Tel +46 8 587 254 00
Fax +46 8 587 255 00
info.europe@elekta.com

Asia Pacific

Tel +852 2891 2208
Fax +852 2575 7133
info.asia@elekta.com

MOSAIQ®

Practice Management

Version 2.5 Features

Overview

Practice Management - Overview	MQ 2.0	MQ 2.1	MQ 2.2	MQ 2.3	MQ 2.4/1	MQ 2.5
User Defined Schedules	✓	✓	✓	✓	✓	✓
Series Scheduling	✓	✓	✓	✓	✓	✓
Appointment Calendars	✓	✓	✓	✓	✓	✓
Scheduling Guidelines	✓	✓	✓	✓	✓	✓
Custom Report Launcher	✓	✓	✓	✓	✓	✓
User Defined Procedure Labels	✓	✓	✓	✓	✓	✓
Multiple Patient Accounts	✓	✓	✓	✓	✓	✓
Batch Edit Procedure Codes	✓	✓	✓	✓	✓	✓
Native PDF File Support	✓	✓	✓	✓	✓	✓
Integrated CodeCheck	✓	✓	✓	✓	✓	✓
Require IDA in demographics	✓	✓	✓	✓	✓	✓
Code Update Enhancements	✓	✓	✓	✓	✓	✓
Require "Responsible" field for QCL	✓	✓	✓	✓	✓	✓
Auto charge enhancements	✓	✓	✓	✓	✓	✓
Weekly Treatment Audit enhancements	✓	✓	✓	✓	✓	✓
Inactive item warning	✓	✓	✓	✓	✓	✓
Require "Account" field in code capture	✓	✓	✓	✓	✓	✓
Wave Edit Enhancements	✓	✓	✓	✓	✓	✓
Kiosk Check-in	✓	✓	✓	✓	✓	✓
ESI Lab Analyzer	✓	✓	✓	✓	✓	✓
ESI Schedule Import	✓	✓	✓	✓	✓	✓
ESI Schedule Export	✓	✓	✓	✓	✓	✓
Enhanced ESI Charge Export	✓	✓	✓	✓	✓	✓
"Mark as Ended" schedule status	✓	✓	✓	✓	✓	✓
Code Set Appointments	✓	✓	✓	✓	✓	✓
UDS Secondary Filter	✓	✓	✓	✓	✓	✓
UDS My Departments View	✓	✓	✓	✓	✓	✓
UDS Quick Appointment Section	✓	✓	✓	✓	✓	✓
UDS Pending Appointment Section	✓	✓	✓	✓	✓	✓

Continued on next page.

Why MOSAIQ?

- Customizable scheduling views to quickly find appointments and optimize resource allocation
- Access patient demographic, clinical and billing information from the schedule to streamline processes
- Intelligent code capture, verification, audit, and export capabilities

For more information please visit

www.elekta.com/MOSAIQ

ELEKTA

Overview Continued

Practice Management - Overview	MQ 2.0	MQ 2.1	MQ 2.2	MQ 2.3	MQ 2.4/1	MQ 2.5
UI Configuration Changes	✓	✓	✓	✓	✓	✓
Additional User Defined Schedule Columns and Options		✓	✓	✓	✓	✓
UDS "Linked" Schedule option		✓	✓	✓	✓	✓
Additional Co-Payment Security		✓	✓	✓	✓	✓
Code Capture History		✓	✓	✓	✓	✓
Enhanced UI Security and Push		✓	✓	✓	✓	✓
ESI Auto Import of Treatment Documents		✓	✓	✓	✓	✓
Charge-Appointment Linking			✓	✓	✓	✓
Staff QCL and Quick Schedule Chart UI			✓	✓	✓	✓
HITECH Meaningful Use (MU) Certified				✓	✓	✓
Login with Windows Authentication				✓	✓	✓
Patient Worklist Enhancements				✓	✓	✓
Add QCLs from the Patient Worklist				✓	✓	✓
Report for QCL Patient Reminders				✓	✓	✓
Patient Demographics Enhancements for MU				✓	✓	✓
Meaningful Use Compliance Reports				✓	✓	✓
Audit Log Enhancements				✓	✓	✓
Batch Edit for Diagnosis Code Library				✓	✓	✓
Diagnosis Code Type Indicator				✓	✓	✓
Demographics Workspace					✓	✓
Authorizations					✓	✓
Exported/Posted Charges Report					✓	✓
Facility Dialog Box					✓	✓
Name Search without Commas					✓	✓
User Defined Schedule Enhancements					✓	✓
Maximum Appointments of Activity Guideline					✓	✓
Biometric Support					✓	✓
Automatic Patient Verification					✓	✓
Quality Check List Enhancements						✓
Episodes of Care						✓
New Patient Registration Window						✓
Demographics Workspace Enhancements						✓
Diagnoses and Problems list pane						✓
Inline Editing						✓
Patient E-mail Column on the User Defined Schedule*						✓
Queuing Across Departments						✓
New Eligibility Verification Work List*						✓
Electronic Eligibility Verification*						✓

*This functionality will be available in the next milestone release of MOSAIQ.

New in MOSAIQ 2.5	Advantage
Quality Check List Enhancements Provides full QCL functionality from the Home and Chart workspace with enhanced fields and filters to help staff find and prioritize tasks.	Improves task management and streamlines QCL maintenance and creation.
Episodes of Care Use "episodes" to record data for individual episodes of care. Use the new Episodes pane in the Demographics workspace to see, add, and change information about the episodes of the global patient.	Define episodes of care and collect episode user defined information.
New Patient Registration Window The enhanced New Patient Registration window has a new look and new fields.	Improved usability.

Demographics Workspace enhancements

New Insurance pane and patient insurance eligibility as well as improvements to the existing workspace.

Diagnoses and Problems List Pane

A new toolbar is added along with improved customization tool

Inline Editing

Use this feature to change the information in the editable field in a pane without going to a different pane or window.

Patient E-mail column on the User Defined Schedule

Use the new *Patient E-mail* column to show the home email address of the patient in the UDS (User Defined Schedule)

New Eligibility Verification Work List

The new Eligibility Work List identifies patients needing eligibility verification based on selected criteria such as upcoming appointments, days since last verification check, and codes captured.

Electronic Eligibility Verification

Electronic Eligibility Verification enables a clinic to easily perform routine batch checks and to update patient insurance information before treatment.

Up to date eligibility status information is accessible to staff throughout the clinic.

Improves user navigation and information access.

Ability to change the layout of the pane for improved information access.

Quickly change an editable field in a pane in a workspace without opening the item or easily go from one item to another and quickly change information in editable fields.

Key patient contact information available during scheduling.

Helps staff verify that patients are covered for service before arrival for an appointment and streamlines eligibility verification tasks.

Reduces time spent verifying patients' eligibility for services.

Accurate insurance eligibility status and coverage information improves reimbursement success.

New in MOSAIQ 2.4/1

Demographics Workspace

The Patient Demographic Data dialog box is changed to the Demographics Workspace with new enhancements including Configurable Workspace.

Authorizations

Authorizations has had multiple enhancements:

- Select the referring physician and referred to physician, and add the issue date.
- Configure the default referring physician for code capture as the Authorization referring physician.

Exported/Posted Charges Report

The Exported/Posted Charges report includes new filter and group options.

Facility Dialog Box

The Facility dialog box has a new field that lets you add a website address.

Name Search without Commas

You can now type a space without a comma between the last and first name to select or find a patient or physician.

Advantage

More efficient accessibility to the information you need, when you need it.

Authorization enhancements include new information and improve flexibility.

Filter charges by location, as well as multiple groupings now available.

Greater functionality in the facility dialog box.

More intuitive and natural for specific regions.

User Defined Schedule (UDS) Enhancements

New features were added to the UDS to help you collect data for Meaningful Use measures. These enhancements let you quickly identify missing information about patients. You can collect the information when the patient checks in or during patient care:

- **New Columns -**
Contact Preference, Diagnosis Stage, Ethnicity, Language, Pharmacy Facility, Race, Deleted By, Deleted Date.
- **Enhanced Columns**
Co-Payments – This column now also shows if a co-payment was skipped
You can configure the UDS so that the new options below appear on the UDS
- **Options Menu**
View Chart, View Observation Orders, View Pharmacy Orders, View Demographics, Restore Appointment.
- **Reports**
You can now configure the Dose Site Summary and Patient Demographics Compliance reports to appear on the Reports sub-menu on the UDS Options menu

Maximum Appointment of Activity Guideline

The new Maximum appointments with this activity rule lets you limit the number of appointments for a certain activity or code group.

Biometric Support

MOSAIQ now uses biometric information for patient identification and verification.

You can set up your workstation to automatically open a patient chart with the fingerprint reader.

Automatic Patient Verification

You can configure your workstations to automatically verify patient before radiation treatment.

- Quickly identify and collect missing Meaningful Use information during patient check-in.
- Open the patient's orders, chart and demographics directly from the highlighted appointment.
- View and restore deleted appointments from any User Defined Schedule.

Restrict or prevent selected activities from occurring at specific times. (i.e., no consultations from 4-6pm on Tuesdays and Thursdays).

You can use the fingerprint reader to enroll fingerprints of a patient in MOSAIQ. The enrolled fingerprints are used for patient identification when they walk-in for an appointment

Speeds up patient treatment workflow with patient information available on the screen before treatment.

www.elekta.com

Human Care Makes the Future Possible

Corporate Head Office:

Elekta AB (publ)
Box 7593, SE-103 93 Stockholm,
Sweden
Tel +46 8 587 254 00
Fax +46 8 587 255 00
info@elekta.com

Regional Sales, Marketing and Service:

North America

Tel +1 770 300 9725
Fax +1 770 448 6338
info.america@elekta.com

Europe, Middle East, Africa, Eastern Europe, Latin America

Tel +46 8 587 254 00
Fax +46 8 587 255 00
info.europe@elekta.com

Asia Pacific


Tel +852 2891 2208
Fax +852 2575 7133
info.asia@elekta.com

ELEKTA

Digital Accelerator

Product Data

1.0 Introduction

This product data describes the high energy, dual modality Elekta digital accelerator for radiation therapy. This accelerator has been designed to deliver therapeutic beams of X-rays and electrons for a wide range of conventional and advanced radiotherapy techniques. These technical specifications apply to the following Elekta linear accelerators:

- Precise Treatment System™
- Elekta Synergy® Platform
- Elekta Infinity™
- Elekta Synergy®
- Elekta Synergy® S*
- Elekta Axesse™ *

** Elekta Synergy S and Elekta Axesse are provided with the Beam Modulator™ radiation head, therefore certain limitations in the energy and fieldsize ranges apply.*

For particular linear accelerator specifications on Beam Modulator-based products consult the corresponding product data:

- Beam Modulator, Functional Description document no. 4513 371 0486
- Elekta Axesse Product Data document no. 4513 371 0603.

This specification does not apply to Elekta Compact™.

1.1 This document replaces the following documents:

- Precise Treatment System Specifications document no. 4513 370 14011
- Elekta Synergy Platform Product Data document no. 4513 371 04271.

1.2 Applicable international standards

The specifications of the linear accelerator declared in this document are based on the recommendations of the International Electrical Commission for the declaration of functional performance characteristics:

- IEC 60976: 1989 incl. Amd 1:2000 Medical Electron Accelerators – Functional Performance Characteristics
- IEC 60977: 1989 Medical Electron Accelerators – Guidelines for functional performance characteristics

The coordinate system convention applied in this document is:

- IEC 61217 (2008) Radiotherapy equipment – Coordinates, movements and scales

Radiation leakage and other safety specifications comply with:

- IEC 60601-2-1 (1998) Medical electrical equipment – Part 2-1: Particular requirements for the safety of electron accelerators in the range of 1 MeV to 50 MeV.

2.0 X-ray beams

The digital accelerator basic configuration includes one X-ray energy beam, which may be selected from tables 1 to 3. Optionally a second and third X-ray energy may be added such that no two energies from the same tables are selected.

2.1 Energy – beam quality

CUSTOMER BENEFIT: Wide choice of energy levels to match your treatment needs

CLINICAL BENEFIT: Ensure best dose distribution by selecting appropriate energy for treatment
Increased flexibility and enhanced efficiency of treatment delivery

Table 1

Nominal energy, MV	4	6
Depth of dose maximum, d_{max} , ± 0.2 cm	1.0	1.5
Percent dose at 10cm depth, D_{10} , ± 1 %	63	67.5
Penetrative quality, d_{80} , cm, ± 0.2 cm	5.9	6.6
Quality index, $TPR_{20,10}$	0.62	0.67
Dose rate variable in steps between (cGy/min.)		
Minimum nominal	25	25
Maximum nominal	230	500

Table 2

Nominal energy, MV	6	8	10
Depth of dose maximum, d_{max} , ± 0.2 cm	1.5	1.9	2.25
Percent dose at 10cm depth, D_{10} , ± 1 %	67.5	71	73
Penetrative quality, d_{80} , cm, ± 0.2 cm	6.6	7.3	7.9
Quality index, $TPR_{20,10}$	0.67	0.71	0.73
Dose rate variable in steps between (cGy/min.)			
Minimum nominal	25	25	25
Maximum nominal	500	600	500

Table 3

Nominal energy, MV	10	15	18	25
Depth of dose maximum, d_{max} , ± 0.2 cm	2.25	2.65	3	3.5
Percent dose at 10cm depth, D_{10} , ± 1 %	73	76.5	78.5	82
Penetrative quality, d_{80} , cm, ± 0.2 cm	7.9	8.9	9.5	10.5
Quality index, $TPR_{20,10}$	0.73	0.76	0.77	0.79
Dose rate variable in steps between (cGy/min.)				
Minimum nominal	25	25	25	25
Maximum nominal	500	600	600	500

Key for tables 1 to 3:

- MV is the nominal energy, it is intended as guidance and for identification of the beam quality
- D_{10} is the percent absorbed dose along the central axis at the depth of 10cm in water in a radiation fieldsize of 10 x 10cm at SSD = 100cm, relative to the maximum dose (Reference BJR Supplement 25)
- d_{max} is the depth of maximum absorbed dose along the central axis in a radiation fieldsize of 10 x 10cm at SSD = 100cm
- d_{80} is the depth of 80% of maximum absorbed dose along the central axis in a radiation fieldsize of 10 x 10cm at SSD = 100cm (Ref. IEC 60976:1989)
- $TPR_{20,10}$ is the ratio of the absorbed doses at depths of 20cm and 10cm in a water phantom, measured with a constant source-detector distance of 100cm and a fieldsize of 10cm x 10cm at the plane of the detector (Ref. IAEA TSR 398, IEC 60976:1989).

2.2 Dose rate

Dose rate in Tables 1 to 3 applies at the depth of the dose maximum in water and 100cm SSD with a 10cm x 10cm fieldsize. The dose rate units of cGy/min are equivalent to MU/min as the accelerator is typically calibrated to provide 1cGy per MU in these experimental conditions.

2.3 Spot size

The spot size of the accelerated electron beam at the X-ray target is less than 2mm as full width half maximum (FWHM), typically 1mm.

2.4 Fieldsize

CLINICAL BENEFIT: Maximum flexibility in treatment area.

- The radiation fieldsize is continuously variable from 0.5cm x 0.5cm to 40cm x 40cm in the plane containing the isocenter. Fields larger than 35cm x 35cm are limited in the corners by a circle of 50cm diameter (defined by the primary collimator).
- Fieldsize is defined as the distance between the 50% dose points on the major axis of the field in the plane of the isocenter at the depth of dose maximum.

2.5 Asymmetrical fields

The standard radiation head is provided with two pairs of independent diaphragms capable of defining fully asymmetrical rectangular fields:

- Y diaphragms (outer) define field edges from 0 to 20cm from the central axis of the beam
- X diaphragms (inner) define field edges from 20cm on one side of the central axis of the beam to 12.5cm on the opposite side (total travel of 32.5cm).

2.6 Uniformity of square X-rays fields

Uniformity of the radiation field is measured in the plane perpendicular to the beam axis at a distance of 100cm from the target (SDD) and at the standard measurement depth of:

- 5cm with an SSD of 95cm for 4MV
- 10cm with an SSD of 90cm for 6MV and above.

The specification values provided apply to the flattened area as recommended in IEC 60976.

- Field flatness – maximum ratio of the maximum absorbed dose to the minimum absorbed dose does not exceed:
 - for field sizes 5cm x 5cm to 30cm x 30cm: $\leq 106\%$
 - for field sizes $>30\text{cm} \times 30\text{cm}$: $\leq 110\%$
- Field symmetry – maximum ratio of absorbed doses at points symmetrically displaced from the axis of the beam and within the flattened area:
 - For field sizes $\geq 5\text{cm} \times 5\text{cm}$: $\leq 103\%$
- Maximum ratio of absorbed dose – for all energies the maximum ratio of absorbed dose in the radiation field to absorbed dose on the radiation beam axis in the plane at the depth of dose maximum:
 - for field sizes 5cm x 5cm to 30cm x 30cm: $\leq 107\%$
 - for field sizes $>30\text{cm} \times 30\text{cm}$: $\leq 109\%$
- Deviation of dose distribution of square X-ray fields with angular positions:
 - maximum variation in the ratio of absorbed dose at any point in the flattened area to the absorbed dose on the radiation beam axis at the standard measurement depth for all angular positions of the gantry and beam limiting system is 2%.

2.7 Penumbra (IEC 60976/ 60977 Clause 7.3)

CLINICAL BENEFIT: Tight penumbra brings advantages for targets close to critical structures
Reduced dose to healthy tissue

Penumbra is defined as the distance between the 20% and 80% dose points, measured in the plane containing the isocenter, at the depth of dose 10cm and SSD = 90cm (5cm and SSD = 95cm for 4MV):

- for field sizes 5cm x 5cm to 15cm x 15cm: $<7\text{mm}$
- for field sizes 15cm x 15cm to 40cm x 40cm: $<8\text{mm}$.

2.8 Leakage and stray radiation

IEC 60601-2-1 (1998), Clause 29: Ionizing Radiation Safety requirements. Patient plane is a plane circular surface of 2 meters radius, centered on and perpendicular to the central axis of the beam and containing the isocenter. The leakage dose measurements are averaged over an area no larger than 100cm² and expressed in percent of the maximum absorbed dose (Gy) at isocenter.

- X-ray leakage absorbed dose (Gy) in the patient plane outside the area of the maximum useful beam, does not exceed 0.1%.
- Neutron leakage absorbed dose (Gy) in the patient plane outside the area of the maximum useful beam, does not exceed 0.02%
- X-ray leakage absorbed dose (Gy) outside the patient plane at a distance of 100cm from the path of the electron beam, does not exceed 0.2%.
- X-ray transmission through the diaphragms. The average X-ray transmission through any diaphragm does not exceed 0.5%.
- Neutron contamination (stray neutron radiation) in X-rays beams. The absorbed dose (Gy) from neutron radiation at the isocenter does not exceed 0.05%.

2.9 Wedged fields

CLINICAL BENEFIT: Increased flexibility to improve conformance to tumor shape

CUSTOMER BENEFIT Efficient treatment delivery times
Removes the physical stress from the operator

A motorized wedge filter is integrated in the radiation head that is automatically inserted in and removed from the radiation beam path:

- wedge angles – continuously variable in the range 0° to 60° (by combining an open field with the 60° wedge field)
- maximum fieldsize – 30cm x 40cm (30cm in the wedged direction).

The automatic wedge eliminates the need for manual selection and insertion of individual wedge filters. The insertion of the wedge filter does not obscure either the light field-indicator or optical distance meter systems. The wedge filter orientation is clearly marked on the radiation head cover and is displayed by the optical field defining light.

3.0 Electron beams

3.1 Energy choice

CLINICAL BENEFIT: Increased flexibility and efficiency in treatment delivery.

The digital accelerator can be configured optionally with up to 9 electron energies in any combination from the options listed in table 4.

3.2 Energies – beam quality

The penetrative qualities given in table 4 are measured along the beam axis at SSD = 95cm and for 10cm x 10cm fieldsize. The parameters listed are (IEC 60976 and BJR-25 definition):

R_{80} – depth of the 80% relative dose

R_{max} – depth of the maximum dose

R_p – practical range of electrons

The nominal electron beam energies in table 4 are calculated as:

$$E(\text{MeV}) = 3 \times \text{dose @ 80\% (i.e. } 3 \times R_{80})$$

The variation of R_{80} does not exceed $\pm 1\text{mm}$ from the values listed in table 4.

Table 4

Nominal energy (MeV)	4	6	8	9	10	12	15	18	20	22
R_{90} , cm	1.35	2.0	2.65	3.0	3.35	4.0	5.0	6.0	6.65	7.35
R_{max} , cm	0.85	1.27	1.7	1.97	2.15	2.6	2.8	3.0	3.1	3.2
R_p , cm	2.1	2.9	3.9	4.5	4.9	5.8	7.3	8.8	9.7	10.8
Dose rate variable in 5 steps (cGy/min.)										
Minimum nominal	25	25	25	25	25	25	25	25	25	25
Maximum nominal	600	600	600	600	600	600	600	600	600	600

CUSTOMER BENEFIT: Enables treatment of tumors when control of depth dose is essential.

3. Dose rate

The electron dose rate is adjustable in five fixed steps from 25cGy/min up to 600cGy/min for all energies. Dose rate is measured with a fieldsize of 10cm x 10cm as defined at SSD = 95cm and at the depth of the dose maximum in a water phantom.

3.4 Fieldsizes

The electron fieldsize is defined with electron applicators at SSD = 95. When attached to the radiation head, the electron applicators do not obscure the light field-indicator or optical distance meter systems:

- standard set of four square electron applicators with fieldsizes of: 6cm x 6cm, 10cm x 10cm, 14cm x 14cm and 20cm x 20cm.
- optional applicators
 - square and rectangular applicators with fieldsizes of: 25cm x 25cm, 20cm x 10cm, 16cm x 8cm and 14cm x 6cm
 - set of four cone applicators with 2, 3, 4 and 5cm diameter
 - an electron arc therapy applicator.
- optional electron beam-shaping mould.

All Elekta electron applicators have removable end frames, which are coded and interlocked for enhanced patient safety. These end frames can be replaced by bespoke end frames, custom made with the electron beam-shaping mould to match the shape of the treatment area.

3.5 Field uniformity (IEC 60976/60977)

(IEC 60976/ 60977 Clause 7) The geometric field edge is defined as the 50% isodose contour measured at the standard measurement depth. The phantom surface is positioned at the end of the electron applicator, SSD = 95cm. Flatness and symmetry are defined within the area confined by a line 1cm inside the 90% isodose contour at the standard measurement depth.

The standard measurement depth for energies of 4 and 6 MeV is the depth of maximum dose. The standard measurement depth for energies of 8 MeV and greater is defined as half the depth of the distal 80% dose point as measured on the central axis of the beam.

- Flatness:
 - The maximum distance between the 90% isodose contour (at the standard measurement depth) and the edge of the geometric field does not exceed
 - on the major axes of the field: $\leq 10\text{mm}$
 - on the diagonal axes of the field: $\leq 20\text{mm}$.

The maximum ratio of absorbed dose anywhere in the radiation field at the standard measurement depth to the absorbed dose on the central axis of the beam at the depth of dose maximum does not exceed 1.03 (103%).

$$\text{Flatness ratio} = \frac{D_{\max}}{D_{\text{axis}}} \quad (\times 100\%)$$

- **Field symmetry:**
the maximum ratio of the maximum to minimum absorbed doses (averaged over $<1\text{cm}^2$ of area) at any two points symmetrical about the central axis of the beam and within the area confined by a line 1cm inside the 90% isodose contour at the standard measurement depth does not exceed 1.03.

$$\text{Symmetry ratio} = \frac{D_{\max}}{D_{\min}} \quad (\times 100\%)$$

X-ray contamination in electron beam (IEC 60601-2-1. 1998, Clause 29.2.1).

Defined as the percentage absorbed dose on the central axis at a depth 10cm beyond the practical range of the electron beam; the X-ray contamination does not exceed:

- for energies 4 to 12MeV inclusive: $\leq 2\%$
- for energies higher than 12MeV: $\leq 4\%$.

3.6 Patient plane leakage (IEC 60601-2-1. 1998, Clause 29.3.1.2 2)

CLINICAL BENEFIT: Reduced dose to healthy tissue.

The average absorbed dose due to leakage radiation in the area between a contour line 4cm outside the periphery of the geometrical radiation field and the boundary of the primary collimator projection at SDD = 95cm does not exceed the 1.5% of the maximum absorbed dose on the beam axis.

3.7 Applicator side leakage (IEC 60601-2-1. 1998, Clause 29.3.1.2)

The leakage absorbed dose measured 2cm outside the surface of the volume containing the body of any electron applicator from its distal end to within 10cm of the enclosure does not exceed 8% of the maximum absorbed dose on the beam axis.

4.0 Dosimetry system

4.1 Construction

CLINICAL BENEFIT: Confidence in providing accurate treatment delivery.

Three independent air-vented ionization chambers aligned in series in the beam path. Two of the chambers (channel-1/primary and channel-2/secondary) provide independent monitoring of the delivered dose. Each dose channel is connected to a separate high precision digital electrometer with automatic leakage compensation and temperature and pressure correction.

One monitor unit (MU) is typically calibrated to 1cGy measured at 100cm SSD in a 10cm x 10cm fieldsize, at a depth of maximum dose.

The third ionization chamber is a multi-electrode system that continuously monitors the radiation beam uniformity. It provides the signals for a highly responsive beam steering and energy servo.

4.2 Termination of treatment

- Normal termination is achieved by coincidence of dose set and dose delivered in channel-1
- Safety (interlock) termination occurs in the event of:
 - channel-1 fails, by channel-2, which is set between 2 and 10 monitor units above the value set in channel-1
 - readings from channel-1 and channel-2 differ by more than 2MU plus 5% at any moment during delivery

- calculated beam-on time has elapsed in the beam timer
- beam non-uniformity measured by the third ionization chamber exceeds 5%
- excess or defect of dose rate is detected.

4.3 Reproducibility

(IEC 60976/ 60977 clause 5.1) Maximum coefficient of variation of the ratio of number of dose monitor units and absorbed dose is $\leq 0.5\%$. This is applicable for any dose rate and to every individual segment of all conformal and intensity modulated delivery techniques.

5.0 Accelerator design features

As well as being designed for clinical flexibility, high performance and excellent accuracy, the advanced technology and digital control used in the Elekta accelerator also provide reduced power consumption and low cost of ownership over the lifetime of the linear accelerator.

5.1 Radio frequency (RF) power source

- The RF power is derived from the compact, high-efficiency MG 6370E Magnetron. The MG 6370E was specifically designed for radiotherapy linear accelerators and features a built-in electromagnetic frequency tuner ensuring precise and responsive AFC.
- This magnetron is able to operate under a wide range of temperature and humidity conditions. It is completely powered off when the linear accelerator is not in operation and will be in operating temperature in less than 10min after powering on.
- The magnetron is mounted on the gantry and rotates with it giving very convenient access for rapid maintenance allowing replacements to be completed in just 30mins. It comes with a full non-pro-rata two-year warranty.

5.2 Accelerating waveguide

Elekta's unique accelerating traveling waveguide is optimized for maximum electron beam throughput and efficient use of RF power. In the traveling waveguide there is no need for complex energy switches subject to mechanical failure and the design provides a unique wide and flexible energy choice. Elekta offers a 20-year warranty with the confidence of knowing that it has never been necessary to replace one of our waveguides in clinical use.

5.3 Electron gun

The electron gun used in the accelerator is of diode type. It is a cost effective and easily replaceable component. The electron gun and magnetron are always pulsed simultaneously avoiding dark current leakage problems and the control of the dose rate is achieved by adjusting pulse repetition frequency.

5.4 Focusing coils

A dual system of focusing solenoids ensures high electron beam transmission through the accelerator structure resulting in reduced stray-radiation and efficient use of RF power.

5.5 Beam steering system

Two pairs of orthogonal steering coils ensure electron beam alignment on the target in all modes, as well as under any gantry orientation. In conjunction with a dedicated multi-electrode ion chamber and digitally controlled servo electronics the steering system minimizes beam symmetry changes under all conditions.

5.6 Slalom beam-bending system

Electron beam-bending is achieved using a patented slalom beam-bending system. This three stage bending magnet system is truly achromatic in both the radial and transverse planes. It ensures the best beam transport control eliminating the need of mechanical energy slits and an electron spot size of typical 1mm, full width at half maximum (FWHM), at exit.

The inherent low profile of the bending magnet results in the lowest isocentric height in dual modality machines in the industry. This feature significantly eases the everyday work of the radiotherapist when setting-up patients for treatment.

5.7 Slitless beam energy control

Precise and stable beam energy is ensured with a sensitive servo control of the gun thermo-ionic emission from a multi-electrode ionization chamber in the x-ray beam path. This servo system eliminates the need for energy-defining mechanical slits in the electron beam path, which would result in increased stray X-radiation.

- electron beam energy spread on target $\leq 5\%$.

6.0 Mechanical parameters and indicators

6.1 Isocenter

CLINICAL BENEFIT: Superior clearance increases flexibility in clinical set-ups and makes non-coplanar beams easy to achieve.

CUSTOMER BENEFIT: A lower isocenter results in:
reduced physical stress for patients and operators
faster set-up times
more accurate positioning of the treatment fields.

- Target to isocenter distance: $100 \pm 0.2\text{cm}$
- Isocenter height above floor: 124cm
- Clearance under radiation head:
 - distance from surface of radiation head to isocenter: 45cm
 - horizontal distance from gantry to isocenter: 124cm
 - radiation head size: 62cm diameter.
- Mechanical and radiation isocenter accuracy:
 - maximum variation of the mechanical and radiation beam axis with rotation of the radiation head and gantry (as defined in IEC 60976/60977): $\leq 1\text{mm}$ radius
 - maximum variation of the table isocenter rotation axis from the radiation isocenter (as defined in IEC 60976/60977) $\leq 0.75\text{mm}$.

6.2 Position indicators scale conventions

Indication of angular and linear positions can be chosen from the following scaling and range conventions:

- IEC 60601-2-1
- IEC 61217
- Bipolar (Elekta).

6.3 Gantry rotation

- Range: $365^\circ (\pm 182.5^\circ)$
- Accuracy of angle indicators
 - digital scale $\pm 0.5^\circ$
 - mechanical scale $\pm 0.5^\circ$
- Resolution of angle indicators
 - digital scale 0.1°
 - mechanical 1°
- Speed of rotation: 0 to 1rpm continuously adjustable.

6.4 Collimator rotation

- Range: $365^\circ (\pm 182.5^\circ)$

- Accuracy of angle indicators
 - digital scale $\pm 0.5^\circ$
 - mechanical scale $\pm 0.5^\circ$
- Resolution of angle indicators
 - digital scale 0.1°
 - mechanical 1°
- Speed of rotation: 0 to 3 rpm continuously adjustable.
- Local control: collimator angle position can be adjusted directly at the radiation head with conveniently placed local controls.

6.5 Fieldsize

CLINICAL BENEFIT: Flexibility in treatment area
Reduced delivery times

- Range: 0.5cm x 0.5cm to 40cm x 40cm
- Accuracy of digital fieldsize indicators:
 $\pm 1\text{mm}$ or $\pm 1.0\%$ per field edge diaphragm (whichever is greater)
- Resolution: 1mm
- Speed of diaphragms: 0 to 1.5cm/sec continuously variable.
- Local controls: diaphragm positions can be adjusted directly at the radiation head with conveniently placed local controls.

6.6 Light field-indicator (field defining light)

Light field-indicator is provided that defines the position of the radiation field in both X-rays and electron modes. Cross wires projected in the light field indicate the position of beam central axis. The integrated wedge filter orientation is also displayed.

- Crosswire accuracy (walk-out) at isocenter: $\leq 1\text{mm}$
- Light field-indicator edge to radiation field edge coincidence:
 - 3cm x 3cm to 20cm x 20cm: $\leq 1\text{mm}$
 - >20cm x 20cm: $\leq 1\%$ of fieldsize
 - radiation field is defined at SDD = 100cm at the 50% of the central dose at the depth of the dose maximum.

6.7 Optical distance meter

The optical distance meter indicates the distance from target (source) to patient surface on the beam axis.

- Range: 75 to 170cm
- Accuracy:
 - $\pm 1\text{mm}$ at 100cm
 - $\pm 2\text{mm}$ at 75cm and 125cm
- Resolution: 5mm.

6.8 Mechanical front pointer

The mechanical front pointer indicates the distance from target (source) to patient surface on the beam axis.

- Range: 85 to 100cm
- Accuracy: $\pm 1\text{mm}$ at 100cm
- Resolution: 5mm.

6.9 Anti-collision protection

CUSTOMER BENEFIT: Confidence in protection of patient.

The head of the digital accelerator is fitted with a positive action touch guard that protects against a collision between the radiation head and the patient on the table or any other object. If activated the interlock chain will stop and inhibit any movements of the gantry, head and table. Temporary override action is available to allow removing the collision conditions.

7.0 Patient support system

The patient positioning system is designed for modern treatment techniques where a high degree of precision is required. The table guarantees that the high standards of stability and repeatability demanded by intensity modulated, image guided and stereotactic techniques are achieved, delivering fast and reliable patient positioning. The patient support system satisfies the recommendations of IEC 60976/60977 for functional performance characteristics.

7.1 Motion ranges

CLINICAL BENEFIT: Versatility in treatment design including complex noncoplanar techniques
Easy and accurate positioning of the patient

- Vertical
 - motorized control
 - range: 650mm to 1750mm (with standard tabletop)
 - speed: 0.2 to 4.0cm/sec continuously controlled
- Lateral
 - manual and motorized control
 - range: ± 250 mm (500mm)
 - speed: 0.2 to 4.5cm/sec continuously variable
- Longitudinal
 - manual and motorized control
 - range: 1000mm
 - speed: 0.2 to 4.5cm/sec continuously variable
- Isocenter rotation
 - manual and motorized control
 - range: $\pm 95^\circ$
 - speed: 0.3°/sec to 5°/sec, continuously variable
- Column rotation
 - manual control with electromagnetic brake
 - range: 360° with indent at 0°.

7.2 Table position indicators

- Translational and vertical
 - accuracy: ± 1 mm
 - resolution: 1mm
- Rotation
 - accuracy: $\pm 0.5^\circ$
 - resolution: 0.1°.

7.3 Maximum patient load

200Kg.

8.0 Tabletops

8.1 Tabletop options

CLINICAL/CUSTOMER BENEFIT: Ergonomic design gives unrestricted access and accuracy even at extremes of travel.

The following table top options are available with Patient Support System:

- C-arm top – a versatile treatment top suitable for therapy, localization and simulation. Movable C-arms permit posterior and oblique fields without the need to move the patient.
- iBEAM® – carbon fiber tabletop with unique advanced design virtually eliminating image artifacts
- iBEAM® evo – next generation of carbon fiber tabletop. This tabletop has no metallic components apart from the rails.

8.2 Tabletop rigidity

- Maximum difference in table height near the isocenter between 30kg load at retracted position and 135kg at extended position is <5mm.
- Additional information and specifications on the patient support system and tabletops can be found in the Precise Table – Functional Description document no. 4513 370 2980.

9.0 Movement controls and display

All motorized movements can be operated simultaneously and at variable speed. Several scaling conventions are available for the customer to choose from (see previous section: mechanical parameters and indicators – position indicators scale conventions).

9.1 Handheld movement controllers

Two handheld movement controllers (one standard, second optional) inside the treatment room allow selection of fieldsize, rotation of gantry and radiation head and control of all motorized table movements. The handheld controller also controls the distance meter, positioning lasers and room lights as well as paging of treatment/machine information displayed on the monitors inside the treatment room.

9.2 Table control panels

Two control panels on the sides of the patient table include controls for longitudinal and lateral movements, isocentric rotation and vertical height adjustments. In addition most movements can be released for manual operation. Rotation about the table support column is manually controlled.

9.3 Motion and position indication

- Patient setting-up is simplified by flat panel screens in the treatment room (one standard, second optional) which provide a visible display of all the positional values as well as patient and treatment prescription information.
- The information display on treatment room monitors can be customized by the customer according to local clinical practice.

9.4 Assisted set-up (ASU)

CLINICAL BENEFIT: Reduced treatment times by simple and accurate patient set-up.

ASU moves the gantry, collimator, and beam geometric parameters and table isocentric rotation to the positions specified in the field prescription. The user can configure the digital accelerator so that the ASU function operates both remotely from the control room and from the handheld controllers inside the treatment room.

- Precision
 - angular positions $<0.5^\circ$
 - linear positions $<0.5\text{mm}$.

10.0 Accessories

10.1 Coded shadow tray assembly (standard)

The shadow tray assembly for beam-shaping with shielding blocks is supplied as standard with two slots for removable shadow trays. The outermost tray can be uniquely identified by means of cut-outs at the edge. The digital accelerator recognizes up to 110 tray codes that can be verified as part of the prescription.

- Maximum fieldsize: 40cm x 40cm
- Maximum shielding block height: 119mm
- Distance shielding block to isocenter: 32.8cm
- Standard package includes the following perspex trays:
 - one holed tray with star pattern
 - one holed tray with parallel slotted pattern
 - one solid (blank) tray.

10.2 Coded shadow tray assembly short (optional)

CLINICAL BENEFIT: Unique identification ensures patient safety

CUSTOMER BENEFIT: Ability to shield critical structures.

An alternative to the standard shadow tray, 2.5cm shorter, it gives increased tray to isocenter clearance facilitating isocentric treatments.

- Maximum fieldsize: 40cm x 40cm
- Maximum shielding block height: 94mm
- Distance shielding block to isocenter: 35.6cm.

10.3 Other optional accessories for beam blocking

- Clinical beam block set for low energy (optional)
 - clinical beam block system consisting of 14 clinical beam blocks for use with photon energies 6MV or less
- Clinical beam block set for high energy (optional)
 - clinical beam block system consisting of 14 clinical beam blocks for use with photon energies 8MV or greater, or for any dual/triple energy systems
- Port film graticule
 - port film graticule using wire inserts. Fits into upper position of the coded shadow tray, enabling simultaneous fitment of blocking tray for treatment verification
- Wire mesh tray
 - fabricated using one piece of aluminum, this tray is used especially with the clinical beam block set.

10.4 Electron applicators

Standard and optional electron applicators are described in the previous section: electron beams – dose rates.

11.0 Optional beam quality matching

CUSTOMER BENEFIT: Increased workflow efficiency.

11.1 Factory match

The option of matching one or more new Elekta machines to each other and/or to an Elekta machine already installed. Matching will be conducted to the factory archived data for the installed machine, or current factory data for new orders. Available for machines post serial No. 105520.

- Photons
 - beam energy: the central axis dose ratio to Dmax at 10cm depth is within $\pm 1\%$ of the average value of the group of machines to be matched for the same nominal energy. The measurement conditions are SSD = 100cm, 10cm x 10cm fieldsize and gantry and BLD at 0° of angles.
 - beam flatness: the relative dose value at any point in the flattened area of the field is within $\pm 1\%$ of the average relative dose value in the same point for the group of machines to be matched. The point measurements are averaged over an area of $\leq 1\text{cm}^2$ and performed at SSD = 90cm and 10cm depth in a fieldsize of 10cm x 10cm and 30cm x 30cm (95cm and 5cm for 4MV). Flattened area is defined as IEC 60976/60977.
- Electrons
 - beam energy: the distal 80% isodose on the central axis is within $\pm 1\text{mm}$ of the average for the group of electron beams, 95cm SSD, gantry and BLD at 0° for 10cm x 10cm²
 - beam flatness – no electron flatness matching is offered.

11.2 Customer match

As factory data match (including match tolerance and conditions) except that the match will be carried out including beam data supplied by the customer. Customer beam data must be supplied with order. Applicable to all machine ranges, but subject to review and acceptance of supplied data.

Additional information and specifications on Beam Matching options can be found in Flatness, Beam Matching and Analysis Protocols, document ref: PMI B173.

12.0 Optional USA electron beam quality specifications

12.1 Energy

At an SSD of 100cm and for a 14 x 14cm radiation field, the maximum deviation of the actual value from the declared value of penetrative quality (see table 4) shall be $\leq 2\text{mm}$. Measured at the distal 80% dose point.

12.2 Flatness

This is defined as the maximum variation in percentage from the mean electron intensity at 100cm SSD, at the depth of R85/2 within the central 80% of the radial, transverse and diagonal axes for field sizes from 10 x 10cm to 25 x 25cm and for all electron beam energies. The maximum flatness variation is $\pm 5\%$.

12.3 Symmetry

Symmetry is defined as the maximum percentage variation in the average electron intensity to the radial and transverse halves of the electron field intensity to the radial and transverse halves of the electron field at R85/2 for a 10 x 10cm and 25 x 25cm field at 100cm SSD. Variation shall not be more than 2% within the central 80% of the geometrical fieldsize.

13.0 Optional beam-shaping devices

CLINICAL/CUSTOMER BENEFIT: Active leakage reduction provided by automatic positioning of diaphragms.

13.1 MLCi

- The Elekta MLCi is a hardware and control software integrated multileaf collimator. It replaces the upper X diaphragms allowing for the external dimensions and excellent treatment clearance of the linear accelerator to be maintained. Both linear accelerator and MLCi are operated from a unique integrated control system ensuring optimum delivery in the most complex techniques.
- The MLCi uses 80 leaves (40 pairs) with 1.0cm resolution at isocenter each one traveling independently up to 32.5cm (12.5cm beyond the central axis) covering a full 40 x 40cm field. It includes backup diaphragms that move in the same direction as the leaves and conventional asymmetric diaphragms in Y direction.
- The MLCi supports conformal RT, step and shoot IMRT, dynamic MLC IMRT and dynamic arc techniques such as VMAT. The diaphragms can operate in automatic leaf-tracking mode ensuring a minimum X-ray leakage in the patient plane.
- Additional MLCi information and specifications can be found in the MLCi, product data, document no. 4513 370 1641.

13.2 MLCi2

CLINICAL BENEFIT: Reduced dose to healthy tissue
Efficient treatment of multiple targets using interdigitation.

The Elekta MLCi2 offers all the technical and clinical benefits of the MLCi enhanced with up to 5 times lower leaf transmission and leakage and interdigitation capabilities. Additional MLCi2 information and specifications can be found in MLCi2, Functional Description 4513 371 0798.

13.3 Beam Modulator™

CLINICAL/CUSTOMER BENEFIT: Fully integrated in MLC while maintaining maximum beam angle flexibility.

- Beam Modulator is a hardware and control software integrated multileaf collimator. It replaces both X and Y diaphragms allowing for the external dimensions and excellent treatment clearance of the linear accelerator to be maintained. Both linear accelerator and Beam Modulator are operated from a unique integrated control system ensuring optimum delivery in the most complex techniques.
- The Beam Modulator comprises 80 leaves (40 pairs) with fine 0.4cm resolution and a maximum fieldsize of 16cm x 21cm. It offers full leaf interdigitation and over-centre travel capable of closing the leaves beyond the 21cm field.
- The Beam Modulator™ supports conformal RT, step and shoot IMRT, dynamic MLC IMRT and dynamic arc techniques such as VMAT.
- Additional Beam Modulator information and specifications can be found in the Beam Modulator, Functional Description document no. 4513 371 0486.

14.0 Optional radiotherapy delivery techniques

14.1 High dose rate electron (HDRE)

Certain medical conditions require total body skin irradiation with low energy electrons. Several techniques have been developed for this purpose to achieve a homogenous dose over the body of the patient and keep X-ray contamination of the beam as low as possible. These techniques require the patient to be positioned a considerable distance from the treatment machine (3 to 4 meters), where the dose rate is considerably lower than at the isocenter. HDRE allows the electron dose rate to be significantly increased allowing HDRE dose rates of 3000 MU/min at isocenter.

- Energy range 4 MeV to 10 MeV
- X-ray contamination (at 3 to 4 meters from isocenter) <1%
- Symmetry at isocenter is $\pm 3\%$
- Maximum dose at isocenter: 10000 MU.

14.2 Intensity modulated delivery techniques

CUSTOMER BENEFIT: A wide selection of treatment delivery options to offer personalized patient care.

Elekta digital accelerators with an MLC head (any of the MLC listed in this document) and operated under control software Desktop Pro™ 7.01 and upwards can be optionally licensed to provide the following advanced intensity modulated delivery techniques.

- PreciseBEAM™ Degmental
 - This option enables the digital accelerator to deliver IMRT in segmental (Step-and-Shoot) mode. Each prescribed segment is delivered accurately and reliably at the highest dose rates for superior patient throughput. Accurate and stable beam control ensures a dose/MU accuracy of $\leq 1\%$ or 0.1MU, whichever is greater, which is vital during the sequential delivery of low dose IMRT fields.
- PreciseBEAM Dynamic
 - With the same excellent dose and geometry accuracy and functionalities as segmental, this option enables continuous dynamic movement of diaphragms and MLC leaves during treatment delivery. Support for popular techniques such as “sliding windows”.
 - The Automatic Leakage Reduction function ensures that during treatment delivery the four individually controlled diaphragms are automatically moved tracking the position where they maximize protection against leakage through MLC leaves, whilst not infringing on the field defined by the multileaf assembly
- PreciseBEAM Dynamic Arc
 - In this arc therapy the linear accelerator delivers a constant number of MU per degree of movement. During delivery simultaneous gantry rotation and motion of diaphragm and MLC leaves is permitted. Dose rate and gantry speed can change along the arc and are automatically selected by the control system to achieve the prescribed dose/deg. Multiple and continuous arcs in CW and CCW direction can be delivered.
- PreciseBEAM™ VMAT
 - This license enables the Elekta volumetric intensity modulated arc therapy (VMAT) treatment delivery. VMAT is capable of simultaneous dynamic control of MLC, diaphragms, gantry and collimator. It allows continuously variable MU per degree along the arc and, as in dynamic arc, the control system automatically adjusts all linear and angular speeds as well as dose rate. Multiple and continuous arcs in CW and CCW direction can be delivered.

CUSTOMER BENEFIT: Reduced delivery times allow for increased patient throughput

CLINICAL BENEFIT: Reduced dose to critical structures
Improved MU efficiency reduces dose to healthy tissue.

Additional information and specifications on intensity modulated delivery options can be found in the Desktop Pro™ 7.01 Control System – Product Data document no. 4513 371 0681.

15.0 Optional IGRT techniques

15.1 iViewGT™

CUSTOMER BENEFIT: Simple and efficient technique for patient positioning verification using easy-to-use registration tools.

- iViewGT is a megavoltage electronic portal imaging device (EPID) intended to augment the existing recording and verification systems providing proof the radiation field has been correctly positioned with respect to the target area. Image detection is performed by means of a scintillator plate attached to an amorphous silicon detector. Image processing within a PC workstation offers a range of options for superimposed graphics, multiple exposures and selectively enlarged images through high-speed analysis and enhancement. Image acquisition is possible before, during and after treatment.
- Additional information and specifications on iViewGT™ can be found in the iViewGT Product Data 4513 371 0294.

15.2 X-ray volume imaging (XVI)

CLINICAL BENEFIT: Ability to visualize soft tissue structures and precisely target dose to the tumor whilst protecting surrounding healthy tissues and critical structures

- The X-ray Volume Imaging system provides high quality kV images for target localization and on-line patient position correction at the time of treatment. XVI provides the option for 2D kV image acquisition using PlanarView™ or MotionView™ and 3D VolumeView™ acquisition with integrated image manipulation and registration tools. Workflow and clinical accuracy are further enhanced by the remote automatic table movement function.
- For additional information and specifications on XVI see the XVI Product Data document no. 4513 371 0650.

16.0 Site requirements

This list includes only some relevant site requirements in general terms. For detailed site planning information please refer to Site Planning Reference documentation no. 1008403.

16.1 Electrical

Electrical supply for linear accelerator: Peak power 30 kVA, Radiating 18kVA. Three-phase, neutral and earth. Nominal voltage 380 to 420V, Nominal frequency 50 or 60Hz.

16.2 Water cooling

A supply of cooling water is required that can be configured as a one-pass system or a closed loop. If the hospital is not ordering an Elekta water cooler then the client is required to supply the linear accelerator with cooled water to the following specification:

- temperature of water at input to the linear accelerator between 12° and 20°C
- maximum flow – 30 liters/minute
- maximum (absolute) pressure at the input to the linear accelerator should not exceed 4 bar.

Maximum heat input into the hospital water is approximately 12kW, thus temperature gain of hospital water at 30 liters/minute flow is approximately 6°C.

16.3 Lighting

There should be no lighting on the ceiling or walls within 500mm either side of the isocenter.

16.4 Cable ducting

Cable ducts are required to run from the rear of the accelerator to the control room. Ducting should be set into the concrete floor for this. Smaller ducts are required to run from the linear accelerator gantry to the water cooler (if used) and to the Client Interface Terminal.

16.5 Lifting equipment

An I-section girder with a safe working load (SWL) of 2200kg should be mounted on to the concrete ceiling directly above and parallel to the rotation axis of the gantry (end stops must be fitted if girder is open-ended).

16.6 Room safety and radiation protection

It should be noted that before constructing or modifying any treatment room, the design must have the approval of the National Radiological Protection Authority. Interlocks must be provided by the customer to interface the treatment room with the Linear Accelerator. These include emergency off switches, room door switches, radiation warning lights and a time delay switch. Connection to these and other customer interfaces is via an interface PCB. The PCB is provided by Elekta.

www.elekta.com

Human Care Makes the Future Possible

Corporate Head Office:

Elekta AB (publ)
Box 7593, SE-103 93 Stockholm, Sweden
Tel +46 8 587 254 00
Fax +46 8 587 255 00
info@elekta.com

Regional Sales, Marketing and Service:

North America

Tel +1 770 300 9725
Fax +1 770 448 6338
info.america@elekta.com

Europe, Middle East, Africa, Eastern Europe, Latin America

Tel +46 8 587 254 00
Fax +46 8 587 255 00
info.europe@elekta.com

Asia Pacific

Tel +852 2891 2208
Fax +852 2575 7133
info.asia@elekta.com

X-ray Volume Imaging (XVI) Product Data

Elekta was the first company to introduce kV imaging at the time of treatment with XVI. With its cutting-edge capabilities, XVI brings 2D, 3D and 4D anatomical soft tissue visualization. With fast integrated imaging and registration tools, XVI effortlessly facilitates patient set-up correction to provide confidence in dose placement. XVI is a flexible imaging solution that gives you choices to ensure the most suitable image guided treatment is available for your patient.

This product data describes in more detail, the XVI solution and the suite of sophisticated imaging tools available.

1.0 Image Acquisition and Registration

XVI provides the option for 2D kV image acquisition using PlanarView™ or MotionView™, 3D acquisition using VolumeView™ and 4D acquisition using Symmetry™.

1.1 MotionView™ kV Sequence Imaging (2D)

CLINICAL BENEFIT:

- Real-time evaluation of patient motion, optimizing delivery

A pre-programmed sequence of images is acquired via preset driven protocols. It enables the viewing of intrafraction motion while the patient is in the treatment position, and is suitable for 2D anatomical motion studies.

1.2 PlanarView™ kV Single Exposure (2D)

CLINICAL BENEFIT:

- Aids initial patient set-up
- High quality images at very low dose (0.2cGy) delivery

Single kV image, equivalent to radiographic mode, suitable for orthogonal or stereoscopic imaging using surrogate implanted seeds/clips/bony anatomy for patient set-up. Acquisition is preset based. Images can be manipulated within XVI, and transferred to MOSAIQ® Image Management software for additional image registration features and trend analysis.

1.3 VolumeView™ 3D Volume Acquisition Mode

CLINICAL BENEFIT:

- Fast acquisition in the treatment position enables real time modification
- Visualisation of all soft tissue structures at the time of treatment decreases risk of geometric miss
- Ease of structure identification

A 3D volume of image data is reconstructed from a series of 2D projection images acquired as the linear accelerator gantry is rotated. This image data can be used for verification of patient position and target motion. The number of projection images acquired can be varied within the preset functionality, depending on the image quality required and patient imaging dose that is considered appropriate for the anatomical region being imaged. This flexibility is provided by utilizing full or partial gantry rotations, with the opportunity to select a choice of gantry rotation speeds.

The reconstruction engine that is included in the VolumeView option is based on a FeldKamp backprojection algorithm specifically optimized for X-ray volume (cone-beam CT) acquisition.

1.4 Symmetry™ 4D Volume Acquisition Mode

CLINICAL BENEFIT:

- Volumetric visualization of the moving tumor at the time of treatment, in the treatment position
-

A 4D volumetric image is reconstructed using a unique anatomy based sorting algorithm which identifies the respiratory phase using anatomical information available within the 2D projection images. Each 2D projection image is sorted into a phase based bin, and reconstructed into a volumetric data set. Within the Symmetry functionality, XVI displays all phase based volumetric images in sequence allowing easy visualization and registration of the moving tumor.

1.5 Intra-fraction Imaging Mode

CLINICAL BENEFIT:

- Visualization of target motion during treatment enables correctional strategies to be implemented
-

This mode allows simultaneous kV imaging during MV treatment delivery.

2D intra-fraction imaging allows visualization of images during treatment either as static or fluoroscopic images and provides an opportunity to monitor and correct for motion during treatment.

3D VolumeView intrafraction imaging allows 3D kV acquisition during rotational MV treatment delivery, providing post treatment verification of the patient position. This allows for online and offline correctional strategies to be implemented.

4D VolumeView intrafraction imaging provides verification of real tumor motion as treatment delivery is undertaken.

1.6 Preset Driven Acquisition and Reconstruction

For both 3D and 4D volumetric imaging, flexibility of acquisition parameters is assured through the implementation of preset parameters that can be configured by the user within the software. These parameters include generator settings, required gantry sweep, appropriate field-of-view settings and collimated X-ray field. All kV acquisition protocols are initiated via a single button press on the integrated function key pad.

Fast in-line reconstruction can be selected for maximum workflow efficiency, reconstruction takes place during image acquisition so that the 3D or 4D image is available immediately following acquisition. The resolution of the reconstruction matrix used can be configured by the user. Presets are supplied for 0.5mm, 0.75mm, 1mm and 2mm voxels

1.7 Remote Table Movement

Workflow and clinical accuracy are further enhanced by the remote automatic table movement function. This enables repositioning of Precise Table with sub-millimeter accuracy in three dimensions following a VolumeView acquisition and image registration. This automatic move can be performed from the control area function keypad, or within the treatment room using the hand held controller.

1.8 Patient Aperture

CLINICAL BENEFIT:

- No compromise in patient set-up positioning
-

The patient clearance aperture is 90cm for both imaging and treatment delivery. The XVI hardware folds away against the gantry when not required. This enables techniques such as non-coplanar treatment techniques to be achieved with maximum efficiency.

1.9 Field-of-view

CLINICAL BENEFIT:

- Largest axial image length

Field-of-view is determined as the visible reconstructed VolumeView image. The available field-of-view options are as follows:

- small = 27cm x 26cm
- medium = 41cm x 26cm
- large = 50cm x 26cm

Below are typical acquisition parameters for VolumeView image guidance.

CTDIw measurements were made with a phantom length of 40cm to account for scatter from a cone beam

CT acquisition

1.9.1 Head-and-neck

CLINICAL BENEFIT:

- 40 second acquisition and reconstruction for low dose, efficient volumetric corrections

Parameter	Value
kV	100
mA	10
ms	10
Field-of view	27cm (small)
Filter	None
No. of projections	183
Gantry sweep (arc)	200 degrees
Gantry speed (deg/min)	360
Reconstruction	In-line
CTDIw	0.5mGy (Head phantom: length 40cm)

1.9.2 Chest

Parameter	Value
kV	120
mA	20
ms	20
Field-of view	41cm (medium)
Filter	Bow tie filter
No. of projections	660
Gantry sweep (arc)	360 degrees
Gantry speed (deg/min)	180
Reconstruction	In-line
CTDIw	4.9mGy (Body phantom: length 40cm)

1.9.3 Prostate Seed

CLINICAL BENEFIT:

- Facilitates fast, efficient registration while still providing full 3D anatomical information
- 40 second acquisition and reconstruction for low dose, efficient volumetric corrections

Parameter	Value
kV	120
mA	16
ms	16
Field-of-view	27cm (small)
Filter	None
No. of projections	183
Gantry sweep (arc)	200 degrees
Gantry Speed (deg/min)	360
Reconstruction	In-line
CTDIw	1.5mGy (Body phantom: length 40cm)

1.9.4 Pelvis

Parameter	Value
kV	120
mA	40
ms	40
Field-of-view	41cm (medium)
Filter	Bow tie filter
No. of projections	660
Gantry sweep (arc)	360 degrees
Gantry Speed (deg/min)	180
Reconstruction	In-line
CTDIw	19.9mGy (Body phantom: length 40cm)

1.9.5 Symmetry

CLINICAL BENEFIT:

- Advanced 4D information with low acquisition doses

Parameter	Value
kV	120
mA	20
ms	16
Field-of-view	27cm (small)
Filter	None
No. of projections	975
Gantry sweep (arc)	200 degrees
Gantry speed (deg/min)	67*
Reconstruction	In-line
CTDIw	12.2mGy (Body Phantom: length 40cm)

* Gantry speed is slowed to ensure adequate respiratory cycles are included in the reconstructed 4D data set.

1.10 Typical file sizes for reconstructed images

	Reconstructed image	Projection data *
VolumeView™ (half rotation scan, 360°/min)	20MB	100MB
VolumeView (full rotation scan, 180°/min)	35MB	350MB
Symmetry™ scan	41MB	525MB

* Projection images can be deleted following reconstruction

2.0 Registration Workflows

2.1 VolumeView™ Registration

CLINICAL BENEFIT:

- Utilize physician defined anatomy for daily registration
- Workflow oriented registration ensures minimum time to analyze, correct and treat the patient in the optimal position

Specific anatomy for registration can be selected by utilizing a clipbox (cube) volume or a shaped region of interest. The shaped registration region of interest allows structures imported from the treatment planning system to be utilized for generation of the registration volume.

The following optimized registration workflows are available for efficient, intuitive image guidance.

2.1.1 Automated Bone Registration (based on chamfer matching)

2.1.2 Grey Value Registration (automated soft-tissue matching)

2.1.3 Manual Registration

2.1.4 Automated Seed Registration

CLINICAL BENEFIT:

- Fast and efficient image guidance with no compromise on volumetric anatomy

This provides an optimized algorithm which facilitates fast, efficient registration of implanted markers within a volumetric environment. Markers/fiducials which have been validated are as follows:

- Visicoil fiducial marker (0.35mm x 10mm)
- CIVCO Gold Soft Tissue Markers (0.8 & 0.9 x 3mm)
- CIVCO Coupled Markers (1.0 x 3mm)
- CIVCO Carbon Marker (1.0 x 3mm)

2.1.5 Critical Structure Avoidance

CLINICAL BENEFIT:

- Increased confidence in accurate treatment delivery, by providing assurance that changes in internal anatomy have not put adjacent critical structures at increased risk since the original treatment plan

This workflow allows registration of the clipbox and the shaped registration region of interest on two separate areas of anatomy. XVI calculates the relationship of both areas of anatomy to the proposed correction vectors. User defined tolerances input into the software will automatically prompt the user if the target has moved closer to the critical structures due to anatomical changes.

2.2 Symmetry™ Registration

CLINICAL BENEFIT:

- Fast, efficient visualisation and registration of the moving tumour at the time of treatment

Each reconstructed phase of the respiratory cycle is automatically matched to a static 3D reference image. Following registration an optimized display allows efficient review of all registration results. Correction vectors can be automatically calculated to position the moving tumor to an average position, facilitating a symmetrical dose delivery, alternatively an exhale matching position can be selected. The Symmetry workflow can be utilized in combination with critical structure avoidance.

3.0 Product Performance

3.1 3D VolumeView™ Image Quality

CLINICAL BENEFIT:

- The ability to clearly visualize soft tissue at the time of treatment allows margin reduction for high dose delivery to the target while protecting critical structures

The 3D image quality has been determined using the Phantom Laboratories CATPhan 503 Phantom with contrast resolution module, spatial resolution module and uniformity resolution module. Example parameters used for image acquisition are:

Parameter	Value
kV	120
mA	20
ms	20
Field-of-view	27cm (small)
No. of projections	650
Gantry sweep (arc)	360 degrees
Reconstruction	in-line
Slice width	1.5mm

3.1.1 Low Contrast Visibility

The contrast resolution module contains 8 x 1.5mm inserts (LDPE, polystyrene, air (x2), Teflon, Delrin™, acrylic and PMP)

- low contrast resolution (using LDPE and polystyrene) $\leq 3.0\%$ at 1mm slice thickness

3.1.2 Spatial Resolution

The spatial resolution module contains 21 spatial resolution sections measuring from 1 to 21 line pairs per cm

- spatial resolution ≥ 10 lp/cm at 1mm slice thickness

3.1.3 Uniformity

The uniformity module contains a material within 2% (2HU) of water

- uniformity is within $\pm 1.5\%$ across a 15cm diameter region at 1mm slice thickness

3.1.4 Geometric Accuracy

CLINICAL BENEFIT:

- Confidence in delivery of high dose treatments

Geometric accuracy is the accuracy of distances measured in the reconstructed data against the physical dimension. The measurements are one using the CATPhan 503 Phantom:

- axial geometric accuracy** $\leq 1\text{mm}$
- sagittal geometric accuracy** $\leq 1\text{mm}$

** dependent on operator.

3.1.5 Hounsfield Unit Accuracy

Hounsfield Unit accuracy: ± 40 HU for a 20cm diameter (air and water equivalent) phantom

3.2 3D System Accuracy

CLINICAL BENEFIT:

- Confidence in delivery of high dose treatments

	RMS error	Max. error
Alignment of MV radiation field center to kV isocenter	0.7mm	1mm
Image registration (on XVI)* Auto-registration bony anatomy	0.5mm	1mm
Automatic table correction	0.5mm	1mm
Total clinical accuracy** (image reconstruction, bony anatomy registration, automatic table correction)	1mm	2mm

* subject to planning CT quality and clinical site

** total clinical accuracy is the accuracy has been determined using a geometric phantom

3.3 4D System Accuracy

CLINICAL BENEFIT: Confidence in delivery of high dose treatments

	Max error
Alignment of Symmetry calculated average position to kV isocenter	1mm
Accuracy of Symmetry calculated amplitude	1mm

3.4 2D PlanarView™ Image Quality

The 2D image quality has been determined using the Leeds TOR 18FG phantom with a 1mm copper plate placed on top of the phantom. The parameters used for image acquisition are:

Parameter	Value
kV	120
mA	10
ms	25
Frames	15
Total mAs	3.75

3.4.1 Low Contrast Visibility

The Leeds test object contains contrast objects from 16% to 0.9%

- low contrast resolution is $\leq 2.7\%$

3.4.2 Spatial Resolution

The Leeds test object contains spatial resolution blocks from 0.5 to 5.0 line pairs per mm

- spatial resolution is ≥ 1.4 lp/mm

3.5 2D System Accuracy

The 2D system accuracy is determined by the difference between the XVI show center position on acquired 2D images on the MV radiation isocenter

- 2D system accuracy < 1 mm

4.0 kV Imaging System

Peak power	40kw, iso spec.		
Radiographic kv range	70kVp – 150kVp		
Max mAs	500mAs		
Voltage ripple	typical $< 1\%$ @ 100kvp	certified	UL, FDA, CSA, CE
Rise time	typical 1ms		
Power input	3 phase		
Weight approx.	75kg		
External dimensions approx.	(w) 545mm (l) 360mm (h) 642mm		

4.1 X-ray Tube

- Fan cooled
- X-ray tube housing assembly: total heat storage capacity: 1200kHU
- Cooling rates for anode and housing (HU/min.) 705HU/s
- X-ray tube and housing filtration 2.6mm Al and 0.1mm Cu

4.1.1 Duty Factor

- Nominally two VolumeView™ scans in 15 mins

4.2 Amorphous Silicon Detector

- 41cm x 41cm
- Nominally 5.5fps
- Image matrix 1024 x 1024 x 16 bits

5.0 DICOM

CLINICAL BENEFIT:

- Easy transfer of information
-

- 5.1 DICOM CT Import
- 5.2 DICOM RT Plan Import
- 5.3 DICOM RT Structure Import
- 5.4 DICOM RT Image Import
- 5.5 Optional DICOM CT Export
- 5.6 Optional Automatic DICOM CT Export (to a Configured Destination)
- 5.7 DICOM 4D CT Export
- 5.8 DICOM RT Image Export

See www.elekta.com for latest DICOM conformance statement

6.0 Calibration Phantoms

6.1 Flexmap Phantom

The single ball bearing phantom and table mount and associated software tools are used for routine QA checks, MV isocenter determination and kV system calibration.

6.2 Optional Phantoms

6.2.1 2D: TOR18FG Leeds X-ray Phantom

The test object contains both a line pair/mm section and low contrast discs.

6.2.2 3D: CATPhan 503 Phantom

Includes contrast resolution module, spatial resolution module and uniformity resolution module.

6.2.3 4D: CIRS Dynamic Thorax Phantom or Quasar Respiratory Motion Phantom

4D phantoms which simulate respiratory motion for analysis of Symmetry functionality.

6.2.4 Quasar Penta-Guide

Cube phantom for daily QA checks of 3D VolumeView registration, kV and MV system coincidence, 2D projection images, laser and field light.

6.2.5 Oval Body Annulus (20cm long)

Oval shaped annulus to simulate the non-uniformity attenuation created by a patient torso shape. Used with the CATPhan for Hounsfield Unit calibration.

6.2.6 Water Calibration Tank

Water tank used for multiple level gain calibration, for enhanced image quality.

7.0 Optional Back-up and Archive

CLINICAL BENEFIT:

- *Storage in MOSAIQ® for distributed review*
-

Tape drive back-up and archive solution – LT03 or LT04 tapes which have 400GB or 800GB native capacity respectively. Alternate archiving solution is available via MOSAIQ® Oncology PACS.

8.0 Optional SYNERGISTIQ™ User Interface

Optional user interface which coordinates MOSAIQ and XVI workflow for IGRT for enhanced efficiency. Allows automatic storage of table correction vectors and CT/RT images within MOSAIQ for offline review, registration and trend analysis.

Note: availability subject to local regulatory clearances. Please contact your local Elekta representative or authorised distributor for details.

www.elekta.com

Human Care Makes the Future Possible

Corporate Head Office:

Elekta AB (publ)
Box 7593, SE-103 93 Stockholm,
Sweden
Tel +46 8 587 254 00
Fax +46 8 587 255 00
info@elekta.com

Regional Sales, Marketing and Service:

North America

Tel +1 770 300 9725
Fax +1 770 448 6338
info.america@elekta.com

**Europe, Middle East, Africa,
Eastern Europe, Latin America**

Tel +46 8 587 254 00
Fax +46 8 587 255 00
info.europe@elekta.com


Asia Pacific

Tel +852 2891 2208
Fax +852 2575 7133
info.asia@elekta.com

Agility™

High Resolution Beam Shaping Product Data

Agility™ from Elekta provides high resolution beam-shaping across a 40 x 40cm fieldsize. The unique low leakage integrated MLC head from the comprehensive Elekta beam-shaping portfolio, builds on Elekta's long experience with optical positioning.

Advanced Clinical Features

- 160 Interdigitating leaves with 5 mm width at isocenter
- Fast leaf speed and high precision
- Integrated dynamic leaf guide removing the need for a split field
- Accurate leaf positioning with Rubicon optical technology*

1.0 Beam Shaping Hardware

Attributes: Mechanical

Interdigitation capable	yes
Number of leaves	160
Nominal leaf width projection at iso-center	5 mm
Maximum field size	40 x 40 cm
Minimum recommended field size	0.5 x 0.5 cm
Maximum distance between leaves on same leaf guide	20 cm
Leaf travel over central axis	15 cm
Leaf nominal height	9 cm
Leaf positioning resolution	0.1 mm
Leaf positioning verification method	Optical
Diaphragm overtravel	12 cm

Leaf Positioning

The robust and reliable Rubicon optical positioning system provides valuable real-time assurance of accurate leaf positioning. Agility's Rubicon optical technology offers advanced real-time leaf monitoring and positioning. Ultraviolet light from an LED source produces infrared fluorescence when it falls on the ruby tips of the multileaf collimator leaves. This infrared fluorescence, detected by an infrared camera, is used to reliably monitor and accurately position the leaves and can be viewed in real-time on the linear accelerator's display screen

Attributes: Dimensions / Weight / Speeds

Head rotation	365°
Head weight	420 kg
Radiation head diameter	81.5 cm*
Head to isocenter clearance	45 cm
Head rotation speed for set-up	12°/s maximum
Head rotation speed for dynamic delivery techniques	6°/s maximum
Leaf speed	Up to 3.5 cm/s Combined with leaf guide up to 6.5 cm/s
Diaphragm speed	Up to 9 cm/s

* Maximum swept diameter

Attributes: Wedge

Integrated wedge size	0-60°
Wedge field size	30 (Y _{IEC}) x 40 (X _{IEC}) cm

Attributes: Physics

X-ray to light coincidence	Maximum distance along the major axes between the light field edge and the radiation field edge for centred fields at normal treatment distance: 5x5 cm to 20x20 cm = 1 mm; 20x20cm to maximum square = 1%
Penumbra (80-20%) for centred fields (at 6 and 10 MV)	< 5.5 mm
Leaf tip penumbra variation for 5x5 cm field over the full travel range (for 6 and 10 MV)	< 1 mm
Leaf position accuracy*	1 mm at isocenter 0.5 mm RMS
Leaf position repeatability	< 0.5 mm
Average transmission through leaves**	<0.375%
Peak transmission through leaves**	<0.5%
X-radiation leakage in patient plane outside primary collimator cone region ***	<0.2% max, <0.1% avg
X-radiation leakage outside patient plane (at 1 m) ****	<0.5%

* Measured using a stripe test. Maximum error quoted as maximum positional error in any leaf pair abutment and root mean square for any leaf pair across all abutments.

** IEC 60601-2-1: 2009, clauses 201.10.1.2.103.2.1a/e and 201.10.1.2.103.2.1b/e for peak and average leakage respectively.

*** IEC 60601-2-1: 2009, clauses 201.10.1.2.103.3a and 201.10.1.2.103.3b for maximum and average leakage respectively.

**** IEC 60601-2-1: 2009, clause 201.10.1.2.104.1a

2.0 Integrated Digital Control System

Integrity™ R3.0 is the seventh generation of integrated digital control system which is the monitoring and controlling foundation of Elekta's digital treatment delivery system. It incorporates the robust medical and aerospace operating system LynxOS. This platform powers the digital linear accelerator offering a system optimized for smooth and efficient delivery of advanced treatment and offering a platform supporting future technologies.

Integrity R3.0 supports a wide variety of techniques from conformal, through IMRT to VMAT for conventional and stereotactic applications, providing clinicians with the freedom to choose the treatment techniques best suited to individual patients.

2.1 Digital Control System Hardware

The Integrity R3.0 software is designed to run on the Agility control system rack cabinet. The control system cabinet (800 mm long x 600 mm wide x 700 mm high) is a small profile cabinet housing dual processors that control both the linear accelerator and the MLC and supports a graphical user interface.

Beam Monitor Unit Display Module (BMDM)

The BMDM module shows the delivered dose in MU. It continues to show the number of monitor units delivered in the last field for a minimum of 20 minutes after power is removed.

User Interface

Elekta Integrity R3.0 user interface is accessed via the control console comprising monitor, keyboard and mouse in the delivery system control area. The graphical user interface supports clinical and service functions and is easy to use and learn. Integrity R3.0 can be configured with up to four monitors, two in the control area and two in the treatment room. The information presented on the monitors is customizable by the user according to local requirements and practice.

Function Key Pad

The Function Key Pad in the control area allows the clinical user to perform Assisted Set-Up (ASU) and to initiate, interrupt or terminate delivery.

2.2 Digital Control System Software

The Integrity software supports two modes of operation which are accessible via a graphical user interface. These modes are designed to meet the needs of the clinical or service user.

2.2.1 Clinical Mode

Provides all the tools required by the clinical user to support and monitor treatment delivery.

Receive External Prescription

Allows the user to deliver radiation using parameters defined by an external record and verify (R&V) system and loaded into the digital accelerator using the Elekta iCom protocol. The patient prescription is stored and managed on the external R&V system.

Receive External Prescription mode also allows the user to copy beams into Service Mode. In Service Mode, copied beams can be delivered as part of any local quality assurance checks or procedures, with no effect on the patient treatment history recorded on the R&V system.

View Item Parts

Allows the performance of selected linac components to be monitored during clinical use.

Assisted Set Up (ASU)

Allows the user to automatically position the gantry, collimator, beam geometric parameters and, optionally the table isocentric rotation to the positions specified in the field prescription. Through the Integrity software, the user can set the configuration of the digital accelerator so that the ASU function operates from the control room and/or the treatment room. It is also possible to disable ASU from the control room when the room doors are open. Table ASU is not possible unless it is licensed and enabled.

CLINICAL BENEFIT: Efficiency in patient set up
Elimination of set up errors
Increased patient throughput

System Administration

Supports the following administration functions

- Manage User Access – used to create and delete users and edit access authorization.
- Customize System Settings – allows the user to define what checks the system will perform when a logged on user attempts to perform a task. It supports user definable tolerance tables to be used in Receive External Prescription.
- Linac Customization – Allows the user to turn linac record printing on or off. Also allows the user to change the valid ranges of the digital accelerator parameters that appear on the user interface. This can be used if, for example, there is a restriction in travel of a certain parameter due to room size.

Beam Timer

Automatic calculation of the backup beam time based on the dose rate. Radiation delivery will terminate when the time is reached.

Support for XVI Image Acquisition

Enables movement of the gantry to the image start position without loading a prescription and gives the option to configure the gantry speed- this feature is enabled through XVI 4.5 or higher.

2.2.2 Service Mode

Provides all the tools required to support routine machine maintenance and calibration, as well as supporting fault finding activities.

Calibration Workflow

Integrity R3.0 introduces new service workflows to support

- Optical system calibration
- Leaf bank height and lateral position set-ups
- Leaf and diaphragm radiation calibration
- X-ray to light field coincidence
- Modules for return to work post-maintenance activities

Intellimax™ Connect

Allows a service user to start or stop a remote access session using Elekta Intellimax.

Service QA tools

- Dynamic log file Viewer – the dynamic log files contain parameter data to assist in reviewing the delivery of dynamic beams. records all parameters for all IMRT, dynamic and VMAT prescriptions delivered.
- DICOM to service mode – direct import of plans from a planning system into service mode
- Copy field via iCom – supports import of clinical prescriptions to service mode via clinical mode.

Configuration Utility

Allows the Service User to configure various aspects of the Control system, for example:

- Perform some standard Windows® configuration functions.
- Configure network printers
- Configure server maintenance activities eg. backup and restore data such as configuration settings and calibration information.
- Configure Intellimax™ settings
- Configure the way parameters are displayed in the user interface
- Configure Scheduled Startup

Diagnostic Utility

Supports the Service User in performing fault finding and maintenance activities, for example:

- Windows® diagnostics and maintenance
- View linac records
- Test network connections
- Cabinet diagnostic tests
- Support for remote support

3.0 Treatment Delivery Options

Agility is driven by Integrity Integrated Control System release 3.0.

3.1 X-ray Beams

Integrity supports up to 3 X-ray beam energies which can be delivered using the following delivery techniques.

Static

Square or irregular shaped beams delivered with a static gantry

Wedge

Supports delivery of wedged fields using an automatic, integrated wedge with angles continuously variable in the range 0° to 60° (by combining an open field with a 60° wedged field). The automatic wedge eliminates the need for manual selection and insertion of individual wedge filters.

CLINICAL BENEFIT: Increased flexibility to improve conformance to tumor shape
Efficient treatment delivery time
Removes manual handling concerns
Wedge insertion automatically controlled by the control system reducing operator errors.

Arc

Gantry rotation during delivery with a fixed field shape and constant gantry speed and dose rate

PreciseBEAM™ Segmental

This optional delivery technique, allows the system to deliver IMRT as a number of sequential segments with different digital accelerator parameters and MLC shapes. The MLC is static for each irradiating segment. Between the delivery of each irradiating segment, the following digital accelerator parameters may move or change as prescribed: Diaphragm position, MLC shape, Wedge position, Radiation energy, Dose rate, Gantry angle and the Collimator angle.

Combining irradiating and move only segments enables the precise delivery of complex IMRT techniques. This includes Step- and- Shoot and Skip-and-Scan arc techniques that use sequential irradiating and move only arcs. Accurate and stable beam control ensures a dose/MU accuracy of $\leq 1\%$ or 0.1MU (whichever is greater), which is vital during the sequential delivery of low dose IMRT fields.

CLINICAL BENEFIT: Improved conformance to target structures
Reduce dose to critical structures

PreciseBEAM™ Dynamic

This licensable option, provides the ability to move the MLC leaves (and diaphragms where applicable) during irradiation, at a specified gantry angle.

CLINICAL BENEFIT: Improved conformance to target structures
Improved delivery speed

PreciseBEAM™ Dynamic arc

This licensable option supports simultaneous movement of the gantry, diaphragms and MLC during irradiation. The beam is delivered with a constant number of MU/degree. Doserate and gantry speed can change along the arc and are automatically selected by the control system to achieve the prescribed MU/Degree. Multiple and continuous arcs in CW and CCW directions can be delivered.

CLINICAL BENEFIT: Improved conformance to target structures
Reduce dose to critical structures and normal tissue
Reduced long and short term side effects.

PreciseBEAM™ VMAT

This optional license provides Volumetric Intensity Modulated Arc Therapy. This technique offers simultaneous dynamic control of the MLC, diaphragms, gantry and collimator. It allows continuously variable MU/degree along the arc and the control system automatically adjusts all linear and angular speeds as well as dose rate to ensure the prescribed MU/degree is achieved. Multiple and continuous arcs in CW and CCW directions can be delivered.

CLINICAL BENEFIT: Enhanced or superior shaping of the dose
Reduction of the dose outside of the tumour volume
Reduced patient plane dose
Fast and efficient delivery
Improved MU efficiency
Reduced patient waiting times
Increased patient throughput

Continuously Variable Dose Rate (CVDR)

This feature enables fine resolution changes in the dose rate for all dynamic delivery techniques, in the range 30cGy to 600cGy/min allowing selection of the ideal value, when delivering either a dynamic or VMAT prescription. This makes the delivery of prescriptions smoother and faster when compared to previous discrete dose rates.

CLINICAL BENEFIT: Smoother, faster delivery
Increased patient throughput
Reducing opportunities for Intra-fraction motion

Interdigitation

Integrity supports interdigitation, enabling planning and delivery of island fields.

CLINICAL BENEFIT: Faster delivery
Improved plan quality
Improved conformance

Hybrid Plan Capability

Allows different treatment techniques to be incorporated in a single delivery.

3.2 Electron Beams

Integrity supports up to 9 electron energies.

Additional Options

Integrity is specifically designed to offer simple upgrades by supporting a range of options that secure investment in this technology.

Remote Automatic Table Movement (RATM)

This licensable option allows the user to perform translational table corrections using values automatically sent from XVI or MOSAIQ following an on-line image guided workflow.

Table moves can be made from either the control or treatment room

CLINICAL BENEFIT: Efficiency in patient set up
Elimination of set up errors
Increased patient throughput

Linac Record

This licensable option allows patient and treatment delivery information to be continuously recorded to a printer. This information can be used as a backup to an R&V system or for billing purposes.

Linac Record to File

This licensable option offers the user the option to send patient and treatment delivery information to a network file rather than to a printer.

Table ASU

This option allows the user to request assisted set up of the patient support system either within the treatment room or from the control area resulting in more efficient clinical workflows.

In-Room Monitor, Keyboard and Mouse

This option provides an additional monitor, keyboard and mouse in the treatment room providing access to all the clinical and service functions of the Integrity R 3.0 control system.

Extended Service

This optional software licence provides additional service tools and functionality which includes

- Scheduled start up
- Display log files
- Edit look up tables
- Create/Edit stored beams
- Full Service graphing
- Calculate Reference Dose
- Table Calibration Wizard

Availability

Please, contact your Elekta representative or an authorized distributor who will advise on regional configuration clearance.

Please contact your local representative for details.

www.elekta.com

Human Care Makes the Future Possible

Corporate Head Office:

Elekta AB (publ)
Box 7593, SE-103 93 Stockholm, Sweden
Tel +46 8 587 254 00
Fax +46 8 587 255 00
info@elekta.com

Regional Sales, Marketing and Service:

North America

Tel +1 770 300 9725
Fax +1 770 448 6338
info.america@elekta.com

**Europe, Middle East, Africa,
Eastern Europe, Latin America**

Tel +46 8 587 254 00
Fax +46 8 587 255 00
info.europe@elekta.com

Asia Pacific

Tel +852 2891 2208
Fax +852 2575 7133
info.asia@elekta.com